

## Transforming Community Equipment Programme

### Product Specification

<b><u>PRODUCT CATEGORY</u></b>	Grab Rails	<b><u>SPECIFICATION NUMBER</u></b>	GR01
<b><u>PRODUCT NAME</u></b>	Drop Down Rail with Leg	<b><u>PRODUCT CODE</u></b>	

#### **FUNCTION:**

- A wall mounted hinged rail with support leg designed to fold up to a wall, locked into position, when not in use and fold down so that the leg rests on the floor for use. The rail is designed to assist the user to stand up or sit down for example onto a toilet pan. The leg must be adjustable in height

#### **MATERIAL:**

- Coated Steel rail and leg
- Where appropriate all materials used must be recyclable

#### **DIMENSIONS:**

- Leg Height – approx 700mm to 900mm.
- Rail Length – approx 750mm
- Base plate – approx 300mm x 125mm

#### **FITTING / FIXTURE REQUIREMENTS:**

- Fixing attachments to be included
- Clear fitting, safe use and cleaning instructions to be included. To include clear instructions regarding fitting on different wall types
- Must comply with Medical Devices Regulations 2002 if applicable
- Where appropriate must be labelled in accordance with CE marking regulations

#### **SERVICING & MAINTENANCE REQUIREMENTS:**

- Not applicable